


Mobilimb™ H3
Hand CPM



**A New
Way
Of
Thinking**



TM Toronto
Medical

Mobilimb™ H3 Hand CPM

- Extremely lightweight
- Portable
- Full composite flexion through 'Spirotary' motion
- Intrinsic Minus and Intrinsic Plus positioning
- Traction treatment
- Several finger-attachment options
- Custom-splint mounting option



Interchangeable springs vary traction force for Traction Treatment.

Several finger-attachment options with easy-to-use Accessory Finger Clips.

The Mobilimb H3 Hand CPM unit demonstrates a new way of thinking. By combining advanced technology with design-simplicity, Toronto Medical has created a hand CPM unit that offers medical professionals more.

Giving you more of what you want.

For instance, more motion and more treatment options through advanced 'Spirotary' motion. This latest innovation in hand CPM enables the H3 to achieve extension, full composite flexion and, for the first time, Intrinsic Minus and Intrinsic Plus positions. All of which means improved range of motion recovery and enhanced tendon gliding.

There are also more treatment application options. Accessory finger clips give you a solution to almost

every finger-attachment situation. We have even included special built-in allowances for applications using custom splints.

For increased treatment flexibility, you can isolate a range of motion anywhere within full range of motion parameters. There is also traction treatment — never before available in an off-the-shelf hand CPM device.

Simple color-coded buttons and dials make adjusting range of motion exceedingly easy. And for quick and easy set up, all intricate adjustments are simplified to the most basic steps. Most people will feel confident applying the H3 after a few simple instructions.

More of what your patients want.

Your patients will appreciate the lightweight, portable design of the H3. Weighing less than 1 lb, the H3 actuator is truly lightweight. Furthermore, with its batteries and electronics conveniently housed in a pocket-size patient controller, it is less obtrusive than other devices.

Of course the H3 is safe for your patients. There is Reverse-on-Load safety circuitry, a conveniently located On/Off button, and a low-battery indicator.

Spirotary Motion Curve



Full extension



Full composite flexion




Intrinsic Plus position

My, how thinking has changed.

In the Mobilimb H3 hand CPM unit less weight, less complication, less nonsense equals more treatment options, more flexibility and

more performance. The Mobilimb H3 Hand CPM unit is more than a new CPM device — it's a new way of thinking.



Color-coded range of motion limiting rings are easy to set and include a safety lock-out mechanism.



Attach the H3 to a custom splint or a cast using the H3 Custom Mounting Kit.