
All the CPM Controller you need.

Mobilimb[™]MFC
Multifunctional Controller





Mobilimb™*MFC*Multifunctional Controller

- Comprehensive control center
- Safety lock-out feature
- Ability to store vital treatment data
- CPM-DATAWORKS Software Package
- Generates a comprehensive Patient-treatment Report

The Mobilimb™ MFC Multifunctional Controller is a comprehensive CPM control center that fits into the palm of your hand. With it you can fully operate the Mobilimb L4 CPM unit.

The MFC allows you to easily and precisely set range of motion, speed and synchronize Muscle Stimulation.

Furthermore, patient safety is enhanced by a highly-accessible Run/Stop button and a safety lock-out feature that prevents settings from being accidently changed.

More than a controller

But there is more to the Mobilimb *MFC* than its function-controlling abilities.

It is also a sophisticated data storage unit that stores all vital patient-treatment data - range of motion, speed, duration of treatment and progress of treatment.

You can then transfer this data to a personal or laptop computer using the Toronto Medical CPM-DATAWORKS Software Package. Once data is transferred it can be viewed on-screen and manipulated to produce a comprehensive Patient-treatment Report.

The Mobilimb *MFC:* The multifunctional controller that puts the control you need plus the versatility you want into the palm of your hand.





