

Mobilimb™ L4
Lower Limb CPM Unit



Setting the standard for Lower Limb CPM



The new Mobilimb™ L4 Lower Limb CPM unit represents a milestone in CPM device technology.

By taking a radically innovative design approach of driving the leg at the knee, instead of the foot, it outperforms conventional CPM devices in several significant areas.

The Mobilimb L4 combines ease of set up, smooth and totally silent operation in a durable, lightweight machine.

Treatment-enhancing features include simultaneous hip, knee and ankle motion from -10° to 135° at the knee, and 4° to 100° at the hip. Furthermore, knee angulation is at a constant speed over the entire range of motion.

Analog controls for range of motion and speed make the L4 easy to operate. Or, for increased operating flexibility an interchangeable programmable remote controller is available.

Other convenient features include a patient control switch and manual fast forward/reverse controls.

Goniometers on the ankle and knee hinges allow for quick and accurate confirmation of angle settings.

And for safety, patented Reverse-On-Load safety circuitry reverses the direction of travel if the patient sufficiently resists motion or if motion is obstructed.

The Mobilimb L4 sets the standard in CPM device technology.

To measure its performance for yourself, call 1 (800) 289-5139 or internationally 1 (416) 420-3303 requesting the Mobilimb L4 information package or to arrange a demonstration.



RANGE OF MOTION FROM -10° to 135° .

TM Toronto
Medical

Helping Patients
Get Better Faster